



# Precellys®裂解试剂盒

不同材质，不同直径的研磨珠，组合超过25种裂解试剂盒可满足您任何样本的需求



产品名称	0.5mL管子	2mL管子	2mL-加厚管子	7mL管子	15mL管子	珠编
Hard tissue grinding <b>MK28</b>	-	P000910-LYSK0-A.0	P000917-LYSK0-A.0	-	-	P000925-LYSK0-A.0
Hard tissue homogenizing <b>CK28</b>	-	P000911-LYSK0-A.0	P000916-LYSK0-A.0	P000935-LYSK0-A.0	P000947-LYSK0-A.0	P000926-LYSK0-A.0
Soft tissue homogenizing <b>CK14</b>	P000933-LYSK0-A.0	P000912-LYSK0-A.0	-	P000940-LYSK0-A.0	-	P000927-LYSK0-A.0
Tough micro-organism lysing <b>VK05</b>	P000934-LYSK0-A.0	P000913-LYSK0-A.0	-	-	P000961-LYSK0-A.0	P000928-LYSK0-A.0
Micro-organism lysing <b>VK01</b>	-	P000914-LYSK0-A.0	-	-	-	P000929-LYSK0-A.0
Soil grinding <b>SK38</b>	-	P000915-LYSK0-A.0	-	P000936-LYSK0-A.0	-	-
Tissue homogenizing CKMix	-	P000918-LYSK0-A.0	-	-	-	-
Bacteria lysing <b>CK01</b>	-	P000919-LYSK0-A.0	-	P000938-LYSK0-A.0	P000948-LYSK0-A.0	P000932-LYSK0-A.0
Micro-organism lysing VKMix	-	-	-	P000937-LYSK0-A.0	-	-
Empty tubes	-	P000943-LYSK0-A.0	-	P000944-LYSK0-A.0	P000946-LYSK0-A.0	-
Tissue grinding <b>CKMix50</b>	-	-	P000922-LYSK0-A.0	P000939-LYSK0-A.0	P000949-LYSK0-A.0	P000930-LYSK0-A.0
Dry hard tissue grinding <b>CK68</b>	-	-	P000923-LYSK0-A.0	-	P000950-LYSK0-A.0	P000931-LYSK0-A.0



样本类型及对应的裂解kit	0.5mL管子	2mL管子	2mL-加厚管子	7mL管子	15mL管子
软组织样本 Starter box P002476-LYSK0-A.0	Soft tissue homogenizing CK14	Hard tissue homogenizing CK28 Soft tissue homogenizing CK14 Tissue homogenizing CKMix			
坚韧组织样本 Starter box P002478-LYSK0-A.0			Hard tissue homogenizing CK28 Hard tissue grinding MK28 Tissue grinding CKMix50 Dry hard tissue grinding CK68		Dry hard tissue grinding CK68
较大组织样本 Starter box P000957-LYSK0-A.0				Soil grinding SK38 Tissue grinding CKMix50 Soft tissue homogenizing CK14	Hard tissue homogenizing CK28
微生物 Starter box P002475-LYSK0-A.0		Tough micro-organism lysing VK05 Micro-organism lysing VK01 Bacteria lysing CK01		Micro-organism lysing VKMix	Bacteria lysing CK01



# 为您的样本选择合适的Precellys 裂解试剂盒

微生物 & 藻类

动物 & 人 - 鱼 & 昆虫

植物 & 土壤

	裂解				匀浆			研磨			
	VK01	CK01	VK05	SK38	CK14	CK28	CKMix	MK28R	CK28R	CK50R	CK68R
藻类			√	√	√	√	√	√	√	√	
细菌 (Gram +/-)	√	√	√	√							
蓝细菌	√	√	√	√							
真菌			√	√							
霉菌			√	√							
浮游植物			√	√							
孢子	√	√	√	√							
病毒	√	√	√	√							
酵母			√	√							
脂肪组织					√	√	√		√		
动脉					√	√	√		√		
骨骼						√	√		√		
脑组织					√	√	√		√		√
软骨					√	√	√		√		√
结肠					√	√	√		√		√
角膜					√	√	√		√		√
耳朵					√	√	√		√		√
粪便				√	√	√	√		√		√
固定过的样本					√	√	√		√		√
头发						√	√		√		√
心脏					√	√	√		√		√
昆虫				√	√	√	√		√		√
肾脏					√	√	√		√		√
幼虫					√	√	√		√		√
肝脏					√	√	√		√		√
肺					√	√	√		√		√
哺乳动物细胞					√	√	√		√		√
肌肉					√	√	√		√		√
皮肤				√	√	√	√		√		√
脾脏					√	√	√		√		√
尾巴					√	√	√		√		√
肿瘤					√	√	√		√		√
蠕虫				√	√	√	√		√		√
叶绿体	√	√	√	√							
棉花							√	√	√	√	√
叶子				√	√	√	√		√		√
苔藓						√	√		√		√
松针						√	√		√		√
植物细胞					√	√	√		√		√
根	√	√	√	√		√	√		√		√
沉积物	√	√	√	√		√	√		√		√
种子				√		√	√	√	√	√	√
土壤				√		√	√		√		√
茎						√	√		√		√
木材						√	√	√	√	√	√



**管子**

- 0.5mL
- 2mL
- 2mL-R
- 7mL
- 15mL


**Prefilled DNase/RNase-free certified kits (per 50 preps)**

**珠子**

(40g to 500g)

- 陶瓷的
- 玻璃的
- 不锈钢的

直径  
from 0.1mm  
to 6.8mm



Bertin中国区总代理: 艾美捷科技  
400-6800-868  
www.amyjet.com